

Plankton Toolbox User's Guide



For Plankton Toolbox version 1.3.2

10 October 2019

Authors: Bengt Karlson and Ann-Turi Skjevik

The Plankton Toolbox Team:

Arnold Andreasson, Marie Johansen, Maria Karlberg, Bengt Karlson and Ann-Turi Skjevik

Illustration on cover page: photo by Ann-Turi Skjevik

Citation

Karlson, B., Skjevik, A-T. 2019. Plankton Toolbox User's Guide version 1.3.2

. Swedish Meteorological and Hydrological Institute, Gothenburg, Sweden, 36 pp.

Plankton toolbox Users Guide

Table of contents

Introduction.....	3
LifeWatch.....	3
Plankton Toolbox is free software.....	3
Uses for Plankton Toolbox	4
Mac, Windows and Linux	5
Getting the latest version of the software and the user's guide	5
Installing the software.....	5
The Plankton toolbox folders	6
Getting help	6
Basic concepts	7
Data sets	7
Taxonomic hierarchy	7
Cell volumes and trophic types	7
The user interface.....	7
Welcome	11
Getting started with Plankton counter	12
Plankton counter	13
Sample info.....	16
Counting methods	17
Count sample.....	18
Sample data	19
Dataset manager	19
Data files – Text(*.txt)	20
Dataset screening	23
Select datasets.....	23
Check structure	24
Check species	24
Check column values	25
Plot parameters	25
Dataset reports.....	26
Dataset analysis.....	26

Select dataset(s)	26
Clean up.....	27
Aggregate/complement data	27
Filter	28
Predefined graphs	28
Generic graphs	29
Statistics	30
Exports.....	30
Managing counting lists.....	30
Managing general species lists.....	31
Important lists	31
Getting started on Mac	32
Technical information for developers	35
Acknowledgements	36
References.....	36

Introduction

Plankton form the base of the food web in most aquatic ecosystems. There is a need to estimate the biomass, abundance and the biodiversity of plankton organisms. Eutrophication, climate change, invasive species and harmful algal blooms are some of the reasons to monitor plankton. Microscope based methods are currently the standard in several monitoring programs including HELCOM-COMBINE, for the Baltic Sea, and OSPAR-JAMP, for North Eastern Atlantic Ocean covering the area between the Azores and the Arctic Ocean. Phyto- and zooplankton samples are collected using e.g. water sampling devices, hoses or nets. Data have been collected for decades and large data sets are available e.g. at international and national data centres. To work with the data in a consistent way may be difficult without the right tools.

The Plankton Toolbox is a free tool for aquatic scientists, and others, working with phyto- and zooplankton data. It is available for MacOS and Windows. Plankton Toolbox makes it relatively easy for non-programmers to work with large data sets on the diversity, abundance, biovolume and carbon content of plankton efficiently. The software is useful for working with datasets emanating from quantitative and qualitative analyses of phytoplankton and zooplankton. Phytoplankton, including harmful algae, are enumerated and identified in numerous ways; see e.g. Karlson et al. (2010). One of the most popular quantitative methods is water sampling, preservation of the sample and subsequent microscope analysis using the sedimentation chamber method (Utermöhl, 1958; Edler and Elbrächter 2010). The method produce data on the biodiversity of plankton. The cell volume of the taxa is also often included to facilitate the calculation of biomass. Plankton toolbox offers a work flow for calculating biovolume of organisms based on Olenina et al. (2006) and also carbon content based on the algorithms by Menden-Deuer and Lessard (2000).


LifeWatch

The development of the Plankton Toolbox is part of the Swedish LifeWatch project funded by the Swedish Research Council Vetenskapsrådet. LifeWatch is also a European Union programme on sharing biodiversity data across Europe.

Plankton Toolbox is free software

Plankton Toolbox is free software and comes with ABSOLUTELY NO WARRANTY.

Uses for Plankton Toolbox



The software Plankton Toolbox has many features, here follows the main ones:

1. The **Plankton counter** provides
 - a. A tool for the microscopist analysing (counting) plankton samples in a consistent way.
 - b. A graphical user interface designed for counting samples using a computer by the microscope. If you prefer pen and paper by the microscope you may find the Plankton counter module useful anyhow.
 - c. A way to record metadata such as sampling data, station name etc.
 - d. A way to store templates with some metadata pre entered, e.g. when samples from different dates from a certain station are analysed
 - e. A way to work with lists of organisms
 - i. The HELCOM-PEG list
 - ii. The Nordic Marine Phytoplankton group list
 - iii. A zooplankton list (ZEN – in development)
 - iv. A list of your own choice
 - v. User defined subsets of lists mentioned above
 - f. Information on traits such as trophic type, cell volume, harmfulness etc.
 - g. Easy calculations of cell abundance and biomass based on concentrated volume, counted area, sedimentation chamber size etc.
 - h. A data format for storing results together with metadata and methods used
 - i. A way to save the results as a report
2. With the **Data set manager** you can :
 - a. Select results from the Plankton counter module
 - b. Import data sets in various formats, e.g. data sets downloaded from data centres.
 - i. Text files
 - ii. Microsoft Excel (.xlsx) files
 - c. Combine different data sets
3. With **Data set screening** you can:
 - a. Carry out some quality control of the data, e.g.
 - i. Screen your data
 - ii. Make plots of the raw data

4. With **Data set reports** you can:
 - a. Select what data you want to include in the report
 - i. Export data in various ways, e.g. in special formats
5. With **Data set analysis** you can
 - a. Clean up your data, e.g.
 - i. to exclude some data in you data set
 1. to select certain species
 2. to select time period
 3. to select station
 4. to select depth interval
 - b. Aggregate/complement your data
 - i. If a sample has been counted at a high level of detail, e.g. at the size group level, you may want to aggregate to a higher taxonomic level, e.g. species, order or class level. Also non-taxonomic plankton groupings may be used.
 - ii. Add zeroes – when combining data from several sampling occasions it is often useful to add zeroes for taxa that have not been observed
 - c. Plot your results
 - d. Calculate statistics

Mac, Windows and Linux

Plankton Toolbox is available for Mac and Windows. A Linux version will be made available upon request. The software has been tested extensively on Windows 7 and to a smaller degree on MacOS 10.13.6, High Sierra. The user's guide provides examples, i.e. screen shots, from Windows 7 and Mac.

Getting the latest version of the software and the user's guide

The software and the user's guide may be downloaded from <http://nordicmicroalgae.org/tools>


Installing the software

Download the zip file. Unpack it and place it in a convenient place on your computer. It should also be possible to run the software on a server or on a virtual machine.

Windows

Place the whole folder with the software and supplied subfolders and files e.g. on the C: drive or the D: drive if you have multiple partitions on your hard drive.

Start the software by double clicking on the black and white icon.



In some computing environments you need to be a local administrator of your computer to be allowed to install or run software. Check with your IT-department if this is needed.

Mac

The instructions below should work for most users.

After unpacking the folder downloaded from <http://nordicmicroalgae.org/tools> you need to do the following:

1. Place the program, i.e. “plankton_toolbox_ver_1_3_0” in the logged in user’s folder, e.g. in Peter, if you are logged in as Peter on your Mac. The folder icon looks like a small house.
2. Also place the subfolders “plankton_toolbox_data and “plankton_toolbox_counter” in the logged in user’s folder, e.g. in Peter, if you are logged in as Peter on your Mac.
3. To start the program double click on “plankton_toolbox_ver_1_3_0”.
4. This will start the window for terminal on the Mac. A few seconds later Plankton Toolbox starts. The terminal window will run in the background. The information shown in the terminal window is not of importance for the user of Plankton Toolbox.
5. Note: the first time you start the software you will need to say ok to running this program on your Mac. You will need to modify security the settings in system settings. First try ctrl-clicking on the software the first time you run it to get a question about running the software directly. See also: <https://support.apple.com/en-us/HT202491>.

For further instructions and screen prints concerning Mac, go to page 32.

The Plankton toolbox folders

plankton_toolbox_data

This folder must not be moved or deleted. Essential parsers, species lists etc. reside in the folder. Read more about this in the section on taxonomic lists near the end of the user’s guide.

plankton_toolbox_counter

This folder must not be moved or deleted. In the config folder, the settings of your counting methods and species lists are saved. In the datasets folder your results are saved automatically.

Getting help

The user’s guide is, at present, the only help system for Plankton Toolbox. The user community is encouraged to use the forum at <http://nordicmicroalgae.org/forum> to post questions and answers and to suggest improvements for the software. Also use the e-mail address info@nordicmicroalgae.org for questions.

Basic concepts

Data sets

Plankton Toolbox treats data as datasets. A data set may contain results from one or several samples. Datasets may be combined, e.g. when working in the *Data manager* and *Dataset analysis*.

Taxonomic hierarchy

One of the features of Plankton Toolbox is the ability to aggregate data to different taxonomic levels, e.g. to class level. This requires a taxonomic tree, i.e. a hierarchy. There are two different taxonomic hierarchies supplied with the package. You may also create your own hierarchy. The hierarchies are user selectable, i.e. the user can use a tree of his or her own choice.

1. The taxonomic hierarchy used in Nordic Microalgae, <http://nordicmicroalgae.org>. This is based on AlgaeBase, <http://algaebase.org> (Guiry and Guiry 2015).
2. The taxonomic hierarchy used in the HELCOM-PEG list.

Cell volumes and trophic types

Another feature of Plankton Toolbox is the ability to work with biovolumes of phytoplankton. A list of cell volumes for phytoplankton taxa from the Baltic Sea region based on Olenina et al. (2006) is supplied with Plankton Toolbox. This list is updated yearly by the HELCOM Phytoplankton Expert Group and is available for download at www.ices.dk. The list also includes information on the trophic type of the organisms, e.g. autotrophic, mixotrophic or heterotrophic. The term not specified (NS) is used for cells that have an unknown trophic type. In addition to the standardised lists support for lists handling synonyms and user defined lists is part of Plankton Toolbox. Calculation of carbon content is part of Plankton Toolbox. The equations used for phytoplankton were developed by Menden-Deuer and Lessard (2000).

The user interface

The user interface consists of one or a few window panes. There is a main window pane and window panes called Extra Tools. You may show and hide window panes as you please. If you have a large computer monitor you may choose to have all open. It is also possible to tear off window panes and place them on the same or on another computer monitor. The Extra Tools may easily be moved around by clicking and dragging them. They may float in front of the toolbox, or placed where one wishes, far right or below the work space.

The Extra Tools

Toolbox logging

In this text file activities are logged and errors reported. Keep the window open when importing new data sets to note problems with species names etc.

Dataset viewer

You can see the original data, the filtered data or choose to not see any data for increased speed.


Data may be exported through the clipboard and pasted into other software such as a text editor or Microsoft Excel.

Data may be exported using the save function as text files or xlsx files.


Graph plotter

In this window pane new plots are shown. Plots may be exported in various formats, e.g. jpg and png.


Here follows some examples of how you may configure the user interface.




Main window only





Main window pane with window pane *Toolbox logging* open




Main window with window pane *Graph plotter* open



Main window pane with *Graph plotter* in a separate window.




Tearing off the *Toolbox logging* window pane to a separate window



Docking the *Toolbox logging* window pane below the main window pane.

Welcome



The Plankton Toolbox welcome page. Text and links will guide you to useful information.


Getting started with Plankton counter

To get started counting a plankton sample you need the following information:


1. Meta data about the sample, where it was collected, when, how it was preserved etc.
 - a. You may want to save the data as a *template sample*
2. A list of taxa that you expect to find in the sample. A useful list for the Baltic Sea area is provided (see managing species lists for details)
3. Information on your counting device e.g.
 - a. Volume of sedimentation chamber or filtered volume, e.g. 20 mL
 - b. The diameter of sedimentation chamber or filtered area, e.g. 26 mm
4. Information on your microscope. You may need to use a stage micrometer (a small ruler) to check the diameter of the field of view in your microscope at a certain magnification. The information may look like this:
 - i. 5x objective xx mm
 - ii. 10x objective yy mm
 - iii. 20x objective zz mm
 - iv. 40x objective åå mm
 - v. 100x objective ää mm
5. Information on sample volume and the volume of preservative added is needed to calculate the dilution of the sample that is a result of adding preservative.
6. The next thing to do is to set up a *method* that suits your work. You may want to set up methods for different magnifications and for different counting styles. Some examples:
 - a. 5x whole chamber
 - b. 10x whole chamber
 - c. 20x transect counting (=counting diameters), 20x field of view
 - d. 40x transect counting (=counting diameters)), 20x field of view
 - e. 100x field of view counting (useful for counting autotrophic picoplankton in the fluorescence microscope)
7. Start counting!

Plankton counter


The plankton counter is a tool for counting zooplankton or phytoplankton samples. The module enables you to record your results while analysing.




Choose one of the existing Plankton counter datasets or create a new one by clicking **New dataset...** and naming the dataset. In the dataset several samples may be added (new samples are added by clicking **New sample...**), so if you have a program for a whole year, you may add all stations, dates and depths analysed.




By clicking **Export/import samples (.xlsx)...** you can import samples that are already created that may be stored somewhere else, on a server or elsewhere. Likewise, you may export samples from your hard disk if you want them to be available from other computers or to colleagues.




Delete incorrect or unwanted samples or datasets:





To enter an existing or new count, mark the actual sample name and click **Edit/count sample...** or double click the actual sample name.



By clicking **Export/import toolbox data (.zip)...** you can backup your work or import earlier backups.

Sample info

Start by unclicking the “locked for editing” box, it is just a security for not altering existing work without intending to.

Plankton counter

Dataset: **2019_National_monitoring** Sample: **Anholt E_2019-08-18_0-10m**

☒ Locked for editing

Sample info | Counting methods | Count sample | Sample data

Sampling

Sample name: Anholt E_2019-08-18_0-10m Change name...

Sample id: Anholt E_2019-08-18_0-10m

Sampling: Year: yyyy Month: mm Day: dd Sampling time (UTC): hh:mm

ICES metadata: Year: Country code: Platform code: Series:

Sampling laboratory: Orderer:

Project:

Method documentation: Method ref.: HC-C-C6

Station name: ANHOLT E Station code:

Latitude, degree: 56 minute: 39.99 Longitude, degree: 12 minute: 7.0

Latitude, decimal: 56.66660 Longitude, decimal: 12.11670

Sampler type code: HOS (Hose)

Sample min depth (m): Sample max depth (m): 10 Water depth (m): 55

Sampled volume (L): 6

Net sampling

Net type code:

Sampler area (m2): Mesh size (µm): Wire angle (deg): Tow length (m):

Analysis

Analytical laboratory: SMHI

Analysis: Year: 2019 Month: 10 Day: 03 Today

Analysed by: Marie Johansen

Comments: Chamber 25-8

Clear Copy from sample...

Close plankton counter

Fill in all relevant information about this sample. If you have a template already, just click

Copy from sample...

and find the template. The information filled in is automatically saved and you can proceed to Counting methods. Note that sampled volume is the volume actually sampled with the hose. This volume is not used when calculating, it is only sampling information.

Counting methods

Plankton counter

Dataset: 2019_National monitoring Sample: Anholt E_2019-08-18_0-10m

☐ Locked for editing

Sample info Counting methods Count sample Sample data

Counting methods

Default method setup: Quantitative_inc_biov Zeiss Axiovert 200

Copy values from selected setup Reset to used values for this sample

Methods: **Method values:** **Method steps:** **Method step values:**

A - Utermohli
B - Filtered sa
C - Filtered sa

Sampled volume (mL): 125
Preservative: CLU (Acid Lugol's solution)
Preservative volume (mL): 0.6
Counted volume (mL): 25
Chamber/filter diameter (mm): 26
Method type: Quantitative

A1 - Utermoh
A2 - Utermoh
A3 - Utermoh
A4 - Utermoh

Magnification: 100
Microscope: Zeiss Axiovert 200
Count area type: Chamber/filter
Diameter of view (mm):
Transect/rectangle width (mm):
Transect/rectangle length (mm):
Coefficient for one area: 40 ☐ Calc. by user
Default counting species list: Baltic phytoplankton list
☒ View sizeclass info

Manage counting methods:

Add method/method step... Delete method step(s)... Save as default method setup... Delete default method setup(s)...

Close plankton counter

Fill in/choose the relevant information. Choose “Default method setup”, default lists are found/saved in the folder named counting_methods, note that you need the correct settings for your specific

microscope. Click **Copy values from selected setup** to get the values from the chosen default method setup. If you are in the middle of a count and accidentally change some values, click

Reset to used values for this sample and your values will be restored.

Choose species list for the count (Default counting species list), mark the “view sizeclass info” if this is an analysis including biovolumes. All cells may be filled in with information of your own choice, also the drop down menu.

Count sample

Plankton counter

Dataset: 2019_National_monitoring Sample: Anholt_E_2019-08-18_0-10m

☐ Locked for editing

Sample info Counting methods Count sample Sample data

Method steps

Method step: A1 - Utermohl - 100x Next step

Count area type: Chamber/filter

Count area number: 1 Add count area Remove count area

Coefficient: 90 Lock taxa...

Species

Scientific name: Guinardia flaccida

Full name: Guinardia flaccida

Sp./spp.: Cf.: Size class: 4

Counting

counted: 5 -100 -10 -1 +1 +10 +100 Clear

Qualitative: 1 2 3 4 5 Clear

Resume counting

Comments:

Summary

Select summary type: Counted per taxa/sizes

☐ Sort on most counted

☒ Current method step only

Total counted: 24

Akashiwo sanguinea [1] : 4 (Size: 40x60)

Guinardia flaccida [4] : 5 (Size: 28-32x1)

Mesodinium rubrum [8] : 2 (Size: 45-55)

Proboscidea alata [7] : 3 (Size: 10x400-500)

Pseudosolenia calcar-avis [4] : 10 (Size: ...)

Species lists

Select counting species list: Baltic phytoplankton list

Filter, part of name: guin Clear

Scientific name	Size class	Sflag	Cells	Trophic type
685 Akashiwo sanguinea	1		AU	Size: 40x60, Volume: 3366
686 Akashiwo sanguinea	2		AU	Size: 40x80, Volume: 4486
687 Akashiwo sanguinea	3		AU	Size: 45x90, Volume: 6390
688 Akashiwo sanguinea	4		AU	Size: 35x50, Volume: 1602
2009 Guinardia delicatula	1		AU	Size: 8-12x20-40, Volume:
2010 Guinardia delicatula	2		AU	Size: 8-12x40-60, Volume:
2011 Guinardia delicatula	3		AU	Size: 13-17x40-60, Volume:
2012 Guinardia flaccida	1		AU	Size: 13-17x60-80, Volume:
2013 Guinardia flaccida	2		AU	Size: 18-22x60-80, Volume:
2014 Guinardia flaccida	3		AU	Size: 23-27x80-100, Volume:
2015 Guinardia flaccida	4		AU	Size: 28-32x110-130, Volume:
2016 Guinardia flaccida	5		AU	Size: 33-37x70-80, Volume:
2017 Guinardia flaccida	6		AU	Size: 33-37x100-150, Volume:
2018 Guinardia flaccida	7		AI	Size: 33-37x200-300, Volume:

Save as counting species list... View sizeclass info Delete counting species lists...

Close plankton counter

To start counting you will need species lists. A phytoplankton list included cell volumes is already in the folder, but you can add your own lists. Now you can start your count. Select species list. Find the species you want to count in the species list (right column) by filling any part of genus or species name. Count by clicking the space bar or by marking the buttons on the screen. Uncount by clicking backspace. The species you count end up in the left column, where you can mark species already in the existing count. You can also save the left hand list as a template for later counts. If you want to create a counting species list before you start, you can create a .txt-file and save it in the following folder:

Plankton Toolbox -> toolbox_data\plankton_counter\config\counting_species_lists

You can easily add or remove counted transects or field views.

Plankton counter

Dataset: 2019_National_monitoring Sample: Anholt_E_2019-08-18_0-10m

☐ Locked for editing

Sample info Counting methods Count sample Sample data

Method steps

Method step: A2 - Utermohl - 200x Next step

Count area type: Transects

Count area number: 2 Add count area Remove count area

Coefficient: 373 Lock taxa...

Species

Scientific name: Nitzschia longissima

Full name: Nitzschia longissima

Sp./spp.: Cf.: Size class: 3

Counting

counted: 4 -100 -10 -1 +1 +10 +100 Clear

Qualitative: 1 2 3 4 5 Clear

Resume counting

Comments:

Summary

Select summary type: Counted per taxa/sizes

☐ Sort on most counted

☒ Current method step only

Total counted: 36

Nitzschia longissima [3] : 4/

Pseudo-nitzschia [4] : 32/0

Species lists

Select counting species list: Baltic phytoplankton list

Filter, part of name: nit Clear

Scientific name	Size class	Sflag	Cells	Trophic type
2529 Nitzschia longissima	1		AU	
2530 Nitzschia longissima	2		AU	
2531 Nitzschia longissima	3		AU	
2532 Nitzschia longissima	4		AU	
2533 Nitzschia paleacea	1		AU	
2534 Nitzschia sigmaidea	1		AU	
2535 Nitzschia sigmaidea	2		AU	
2536 Nitzschia vitrea	1		AU	
2537 Nitzschia	1		AU	
2538 Nitzschia	2		AU	
2539 Nitzschia	3		AU	
2540 Nitzschia	4		AU	
2541 Nitzschia	5		AU	
2542 Nitzschia	6		AI	

Save as counting species list... View sizeclass info Delete counting species lists...

Close plankton counter

Click [Next step](#) or open the drop down menu to start counting on the next level, e.g 200x or 400x.

Click [Lock taxa...](#) if you have an organism that you want to stop counting after x transects or field views, but you want to continue counting transects or field views with the rest of the sample.

Sample data

Plankton counter ? ×

Dataset: **2019_National_monitoring** Sample: **Anholt E_2019-18_0-10m**

☐ Locked for editing

Sample info Counting methods Count sample **Sample data**

	scientific_full_name	taxon_class	scientific_name	size_class	method_step	count_area_number	locked_at_area	counted_units	counted_units_list	abundance_class	coefficient	abundance_units_l	volume_mm
1	Akashiwo sanguinea	Dinophyceae	Akashiwo sanguinea	1	A1 - Utermohl - 100x	1		4	4		40	160.0	0.005386
2	Guinardia flaccida	Bacillariophyceae	Guinardia flaccida	4	A1 - Utermohl - 100x	1		5	5		40	200.0	0.01696
3	Mesodinium rubrum	Litostomatea	Mesodinium rubrum	8	A1 - Utermohl - 100x	1		2	2		40	80.0	0.005233
4	Nitzschia longissima	Bacillariophyceae	Nitzschia longissima	3	A2 - Utermohl - 200x	2		4	4;0		373	1492.0	0.0008636
5	Proboscia alata	Bacillariophyceae	Proboscia alata	7	A1 - Utermohl - 100x	1		3	3		40	120.0	0.004239
6	Pseudo-nitzschia	Bacillariophyceae	Pseudo-nitzschia	4	A2 - Utermohl - 200x	2		32	32;0		373	11936.0	0.007305
7	Pseudosolenia calcar-avis	Bacillariophyceae	Pseudosolenia calcar-avis	4	A1 - Utermohl - 100x	1		10	10		40	400.0	0.03561


[Export sample \(.xlsx\)...](#) [Save edited changes](#)

[Close plankton counter](#)

On the sample data page you get an overview of your count. If you detect mistakes you may change species or delete entire posts. When you are ready click [Export sample \(.xlsx\)...](#) and save your results wherever suitable.


Dataset manager

The Plankton Toolbox has a dataset manager where you can import data from different sources to manipulate in several ways.



Plankton Toolbox imports data in text (txt) and Microsoft Excel (xlsx) files in different formats in addition to the plankton counter samples. Since data formatting differs depending on the source of the data, importing formats are configurable by the user through parsers. A useful function is that data from different sources, in different formats, can be combined by importing multiple files. The data can then be exported as one consistent dataset in txt or xlsx format by the user.

Data files – Text(*.txt)



To work with special text- or Excel files with other rows and names than these predefined you may need a custom-made parser which SMHI can help you with. In the text files tab you can for instance manage data downloaded from the Sharkweb, <https://sharkweb.smhi.se/>, at present only available in Swedish.

URVAL PROVER OCH MÄTVÄRDEN

Urval prover och mätvärden Avancerat urval [Hjälp](#) ☐ Dölj urval

När:

År: 2016 - 2018

Månader: [Alla](#) [Rensa](#)

- ☒ Januari
- ☒ Februari
- ☒ Mars
- ☒ April
- ☒ Maj
- ☒ Juni
- ☒ Juli
- ☒ Augusti
- ☒ September
- ☒ Oktober
- ☒ November
- ☒ December

Vad:

Datatyp:

Parameter:

Stationsnamn:

Art/taxon-namn:

Visa på kartan:

Visa som tabell/rapport:

Valj alternativ för visning/nerladdning beroende på datatyp och behov:

Tabell:

Rubrikrad:

Visa som diagram:

Diagram:

Geografisk avgränsning ej aktiverad. Kan aktiveras till höger om kartan.

When downloading from the Sharkweb, make sure you choose “internt namn” (internal name), then click “sök” (search).

RESULTAT PROVER OCH MÄTVÄRDEN (Rad 1 till 2000 av 12110 visas) [Nästa sida](#)

2	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
3	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
4	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
5	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
6	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
7	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
8	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
9	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
10	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
11	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
12	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs
13	Phytoplankton	Klar	Leverantör	2018	ANHOLT E	ANHOLT E	NAT Nationell miljöövervakning	Havs

< >

DIVERSE


Ladda ner data Allmänna inställningar

Decimal/fält-avgränsare: Radbrytning: Teckenkodning:


[Villkor för nerladdning](#) ☒ Jag accepterar villkoren

Resultat med fler än 50000 rader levereras osorterat.


The results will appear, accept the downloading terms (jag accepterar villkoren), and click “ladda ner data” (download). Now your data is ready for the data manager.



Click the import dataset(s) to browse for the data file(s) that you want to work with.




You may add files with other formats, for example from the plankton counter:



Dataset screening


When the data has been imported the raw data can be screened in different ways to check for duplicate data, look for unrealistic dates, positions etc. There are also plotting tools to visualize the raw data.

Select datasets



The dataset(s) imported in the data manager can be screened here, mark the one(s) to be screened.


Check structure




Check the data structure and also get a quality check by clicking



Check species




Check column values



Each column may be marked in the drop down menu and checked.

Plot parameters



Mark the parameter(s) you want to plot and click either **Plot parameter values** or

Plot parameters values per date

Dataset reports

Plankton Toolbox - Desktop application

File Extra tools Help

Activities:

- Welcome
- Plankton counter
- Dataset manager
- Dataset screening
- Dataset reports**
- Dataset analysis

Extra tools:

- Dataset viewer (hide)
- Graph plotter (hide)
- Toolbox logging (hide)
- (hide all)

Dataset reports

Select datasets

☒ Import-1: Source: AnholtE_2016-2018.txt

☐ Import-2: Source: 2019_National_monitoring: Anholt E_2019-08-18_0-10m

Clear all Mark all

Create report

Report type: Quantitative (counted): Table format ☒ Aggregate similar rows

Create report

Report preview

	Station name	Station code	Sample date	Analysis date	Min depth (m)	Max depth (m)	Class	Scientific name	Species flag code	Cf	Author
1	ANHOLT E		2016-01-12		0.0	10.0		Ciliophora			
2	ANHOLT E		2016-01-12		0.0	10.0	Cryptophyceae	Cryptomonadales			
3	ANHOLT E		2016-01-12		0.0	10.0		Flagellates			
4	ANHOLT E		2016-01-12		0.0	10.0	Dinophyceae	Gymnodiniales			
5	ANHOLT E		2016-01-12		0.0	10.0	Raphidophyceae	Heterosigma akashiwo			
6	ANHOLT E		2016-01-12		0.0	10.0	Ciliata	Laboea strobila			
7	ANHOLT E		2016-01-12		0.0	10.0	Katablepharidophyceae	Leucocryptos marina			

Save report

Copy to clipboard File format: Excel file (*.xlsx) Save as...

A handy way to format a dataset into a data report for delivery, for instance for an outsourcer. Click the Aggregate similar rows box if the taxa are not needed to be presented in the different size classes, there will be an aggregation on a taxon level.

Dataset analysis

Select dataset(s)

Plankton Toolbox - Desktop application

File Extra tools Help

Activities:

- Welcome
- Plankton counter
- Dataset manager
- Dataset screening
- Dataset reports
- Dataset analysis**

Extra tools:

- Dataset viewer (hide)
- Graph plotter (hide)
- Toolbox logging (hide)
- (hide all)

Dataset analysis

Select dataset(s) Clean up Aggregate/complement data Filter Predefined graphs Generic graphs Statistics Exports

☒ Import-1: Source: AnholtE_2016-2018.txt

☐ Import-2: Source: 2019_National_monitoring: Anholt E_2019-08-18_0-10m

Clear analysis data Load marked dataset(s) for analysis

Analysis data, filtered data, statistical data and export data

View: Analysis data Number of rows: 12110


	visit_year	sample_date	visit_month	station_name	sample_latitude_dd	sample_longitude_dd	water_depth_m	sample_id	sample_min_depth_m	sample_max_depth_m	scient
1	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Attheya septen
2	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Attheya septen
3	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Attheya septen
4	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Attheya septen
5	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Bacillariophyce
6	2018	2018-12-08	12	ANHOLT E	56.6666	12.1167	55.0	Anholt E_2018-12-08_0-10m	0.0	10.0	Bacillariophyce

Export data

Copy to clipboard File format: Tab delimited text file (*.txt) Save...


To analyse data, select the wanted dataset(s), click Load marked dataset(s) for analysis and the data will appear in the lower window. Observe that the data at any time can be saved in two different ways, either by clicking **Copy to clipboard** or **Save...**, remember to choose format.

Clean up



In the “clean up” step you can remove data you do not need, like certain dates, stations, depths, trophic types etc.

Aggregate/complement data




At this step data can be aggregated to different taxonomic levels, e.g. from species level to class level. There is a function for adding zeros for organisms that are found in some but not all samples. The software looks through all the sampling locations, dates and depths, creates a complete list of taxa observed in all the samples in the dataset, and adds zeros in abundance for taxa that were not observed in certain samples. It is possible to choose whether to include trophic types and development stage/sex (zooplankton).

Filter

Here you can filter your data further, for example to exclude certain species or months.


Predefined graphs

A few graphs are available for quality control and such. With the Extra tools in the left margin you can view/hide the dataset, graphs or loggings. The current aim of the plotting tools is to give the user the ability to produce fairly simple graphs. For publication quality output it is often necessary to export the data and to use some graphical software package.



A time series is one of the predefined graphs.

Generic graphs



This tab offers more options in generating graphs.

Statistics

Plankton Toolbox - Desktop application

File Extra tools Help

Activities:

- Welcome
- Plankton counter
- Dataset manager
- Dataset screening
- Dataset reports
- Dataset analysis**

Extra tools:

- Dataset viewer (hide)
- Graph plotter (hide)
- Toolbox logging (hide)
- (hide all)

Dataset analysis

Select dataset(s) Clean up Aggregate/complement data Filter Predefined graphs Generic graphs **Statistics** Exports

Parameter: Abundance (ind/l or 100 um pieces/l)

Split by:

- ☒ Year
- ☐ Season
- ☐ Month
- ☐ Station
- ☒ Sampling event
- ☒ Depth
- ☐ Scientific name

View data Calculate statistics Plot graph

Analysis data, filtered data, statistical data and export data

View: Statistical data Number of rows: 53

	Parameter	Year	Sampling event	Depth	Mean	Median	Std. dev.	Min	Max	Counted values
1	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-01-12	0.0-10.0	10142.391304347826	1145.0	17736.075071273164	50.0	69069.0	23
2	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-02-19	0.0-10.0	11314.472222222223	1030.5	21216.89864168302	50.0	82509.0	36
3	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-02-21	0.0-10.0	6155.613636363636	1759.0	9361.57677864906	49.0	37972.0	44
4	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-03-18	0.0-10.0	11788.825	2831.0	18732.40024114302	49.0	74132.0	40
5	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-03-19	0.0-10.0	18229.204545454544	4260.0	33381.73245720666	100.0	144840.0	44
6	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-04-21	0.0-10.0	27974.14285714286	5133.0	63776.47213069044	50.0	282315.0	35

Export data

Copy to clipboard File format: Tab delimited text file (*.txt) Save...

Simple statistics may be performed and plotted.

Exports

Plankton Toolbox - Desktop application

File Extra tools Help

Activities:

- Welcome
- Plankton counter
- Dataset manager
- Dataset screening
- Dataset reports
- Dataset analysis**

Extra tools:

- Dataset viewer (hide)
- Graph plotter (hide)
- Toolbox logging (hide)
- (hide all)

Dataset analysis

Select dataset(s) Clean up Aggregate/complement data Filter Predefined graphs Generic graphs Statistics **Exports**

Parameters (for PRIMER):

- ☒ # counted (ind/analysed sample fraction)
- ☒ Abundance (ind/l or 100 um pieces/l)
- ☒ Biovolume concentration (mm³/l)
- ☒ Carbon concentration (ugC/l)

PRIMER

Zooplankton: Abundance m2 and m3, length median and mean

Clear all Mark all

Analysis data, filtered data, statistical data and export data

View: Statistical data Number of rows: 53

	Parameter	Year	Sampling event	Depth	Mean	Median	Std. dev.	Min	Max	Counted values
1	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-01-12	0.0-10.0	10142.391304347826	1145.0	17736.075071273164	50.0	69069.0	23
2	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-02-19	0.0-10.0	11314.472222222223	1030.5	21216.89864168302	50.0	82509.0	36
3	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-02-21	0.0-10.0	6155.613636363636	1759.0	9361.57677864906	49.0	37972.0	44
4	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-03-18	0.0-10.0	11788.825	2831.0	18732.40024114302	49.0	74132.0	40
5	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-03-19	0.0-10.0	18229.204545454544	4260.0	33381.73245720666	100.0	144840.0	44
6	Abundance (ind/l or 100 um pieces/l)	2016	ANHOLT E 2016-04-21	0.0-10.0	27974.14285714286	5133.0	63776.47213069044	50.0	282315.0	35

Export data

Copy to clipboard File format: Tab delimited text file (*.txt) Save...

Arranges the data in for example PRIMER format, a statistical tool specialized for biological data.

At most steps in the data processing data can be copied to the clipboard and pasted into other software running on the computer used. Data may also be exported in txt or xlsx formats for further analyses or plotting using other software.

Managing counting lists

When you count a sample you may click on *Save counting species list as* and save a list of the taxa and size classes you may want to re-use in the future. The counting list will show up in the drop down window to the left in the counting window pane. Keep in mind that the counting lists for counting samples always are subsets of the larger species lists. The subsets are found in the folder:

Plankton Toolbox -> toolbox_data\plankton_counter\config\counting_species_lists

You may change these lists as you please but make sure that you only include taxa and size groups that are found in the general species lists (see separate section on this topic).

Managing general species lists

Plankton Toolbox provides species lists based on the HELCOM-PEG groups lists. You may also use other lists or create your own. Keep in mind that you need to include all organisms found in your sample or your data set in the species lists. Organisms missing in the taxonomic hierarchy will not be included when data is aggregated to different taxonomic levels.

The files with the species list and lists on taxonomic hierarchy are found in the folder Plankton Toolbox -> toolbox_data\species

The lists are imported into Plankton Toolbox when the software is started.

Important lists

1. A list of taxa that you expect to find in the samples you analyse or in the dataset you work with
 - a. The HELCOM-PEG 2019 list is provided. It includes a list of taxa, trophic type, cell shape, cell volume, carbon content per cell etc.
 - b. An amendment to the list focusing on organisms in the Kattegat-Skagerrak may also be included.
 - c. You may add taxa in the file taxa_extra.xlsx file. Always include information on taxonomic hierarchy, e.g. rank and the name of the taxon higher in the taxonomic hierarchy.
2. To be able to aggregate data to different taxonomic levels a taxonomic hierarchy is needed
 - a. In the HELCOM-PEG 2019 list the following taxonomic levels are provided:
 - i. Division
 - ii. Class
 - iii. Order
 - iv. Scientific name (most often species)
 - b. In the Nordic Microalgae list (it will be added at a later stage) the following taxonomic levels are provided:
 - i. Biota
 - ii. Phylum
 - iii. Class
 - iv. Order
 - v. Genus
 - vi. Species
 - vii. Subspecies or Variety or Forma
 - c. In the list taxa_extra.xlsx you may add new taxa with information on taxonomic hierarchy
3. A list defining synonyms and translations from misspelled to correct names is also provided. The file name is translate_worms.xlsx.
4. A list defining harmful species is provided, called harmful_hab. Harmful_extra defines the higher levels, genus or order which are missing in the harmful_hab list.

Keep in mind that the counting lists always are subsets of the larger general species lists. The subsets are found in the folder:

Plankton Toolbox -> toolbox_data\plankton_toolbox_data\species.

Getting started on Mac

Download the zip archive with the software and accompanying files from:

<http://nordicmicroalgae.org/tools>


Unpack the archive, e.g.

Mac


The instructions below should work for most users.

After unzipping the folder downloaded from <http://nordicmicroalgae.org/tools> you need to do the following:

1. Place the program, i.e. "plankton_toolbox_ver_1_3_0" in the logged in user's folder, e.g. in Peter, if you are logged in as Peter on your Mac. The folder icon looks like a small house.
2. Also place the subfolders "plankton_toolbox_data and "plankton_toolbox_counter" in the logged in user's folder, e.g. in Peter, if you are logged in as Peter on your Mac.
3. To start the program double click on "plankton_toolbox_ver_1_3_0".
4. This will start the window for terminal on the Mac. A few seconds later Plankton Toolbox starts. The terminal window will run in the background. The information shown in the terminal window is not of importance for the user of Plankton Toolbox.




You may choose to run Plankton Toolbox in full screen mode. Just click on the green button in the upper left part of the window.



Plankton Toolbox running in full screen mode.

```
Bengt — plankton_toolbox_start — plankton_toolbox_start * plankton_toolbox_start — 80x24
Last login: Mon Dec 7 20:04:48 on ttys000
/Users/Bengt/PlanktonToolbox_ver_1_1_1_Mac_2015-11-07/plankton_toolbox_start ; e
xit;
Bengts-MacBook-Air:~ Bengt$ /Users/Bengt/PlanktonToolbox_ver_1_1_1_Mac_2015-11-0
7/plankton_toolbox_start ; exit;
QPixmap::scaled: QPixmap is a null pixmap
QPixmap::scaled: QPixmap is a null pixmap
QPixmap::scaled: QPixmap is a null pixmap
```

The Window for the Mac Terminal will be running in the background. The information shown in the Terminal window is not of relevance for the user of Plankton Toolbox.




```
Bengt — plankton_toolbox_start — 80x24
Last login: Mon Dec  7 20:04:48 on ttys000
/Users/Bengt/PlanktonToolbox_ver_1_1_1_Mac_2015-11-07/plankton_toolbox_start ; e
xit;
Bengts-MacBook-Air:~ Bengt$ /Users/Bengt/PlanktonToolbox_ver_1_1_1_Mac_2015-11-0
7/plankton_toolbox_start ; exit;
QPixmap::scaled: QPixmap is a null pixmap
QPixmap::scaled: QPixmap is a null pixmap
QPixmap::scaled: QPixmap is a null pixmap
QPixmap::scaled: QPixmap is a null pixmap
logout
Saving session...
...copying shared history...
...saving history...truncating history files...
...completed.

[Processen slutförd]
```

You may want to quit the Mac Terminal after you have quit Plankton Toolbox.

Technical information for developers

The software Plankton Toolbox was developed by Arnold Andreasson using open source software, i.e. Python version 3.x. The code is free as defined by the MIT-license, the Open Source Initiative, <http://opensource.org/licenses/mit-license.php>.



Acknowledgements

The development of Plankton Toolbox was supported by the Swedish Research Council through the Swedish LifeWatch project. The effort by phytoplankton specialists who tested the software and suggested improvements is much appreciated.

References

Edler, L., Elbrächter, M. 2010. The Utermöhl method for quantitative phytoplankton analysis. in: Microscopic and molecular methods for quantitative phytoplankton analysis, Eds. Karlson, B., Cusack, C. and Bresnan, E. Intergovernmental Oceanographic Commission Manual and guides, pp 13-20.

HELCOM 2015 Manual for Marine Monitoring in the COMBINE Programme of HELCOM. 413 pp. World-wide electronic publication, Helsinki Commission, <http://helcom.fi/action-areas/monitoring-and-assessment/manuals-and-guidelines/combine-manual> downloaded on 24 June 2015.

Karlson, B., Godhe, A., Cusack, C. and Bresnan, B. (2010) Overview of methods for quantitative phytoplankton analysis in: Microscopic and molecular methods for quantitative phytoplankton analysis, Eds. Karlson, B., Cusack, C. and Bresnan, E. Intergovernmental Oceanographic Commission Manual and guides, pp5-12.

Karlson, B., Andreasson, A., Johansen, M., Mohlin, M., Skjevik, A-T., Strömberg, P. (2015). Plankton Toolbox – open source software making it easier to work with plankton data, Proceedings of the 16th International Conference on Harmful Algal Blooms, 27 Oct. to 1 November, 2014, Wellington, New Zealand.

Menden-Deuer, S., and Lessard, E. J. 2000. Carbon to volume relationships for dinoflagellates, diatoms, and other protist plankton. *Limnology and Oceanography*, 45: 569–579.

Olenina, I., Hajdu, S., Edler, L., Andersson, A., Wasmund, N., Busch, S., Göbel, J., Gromisz, S., Huseby, S., Huttunen, M., Jaanus, A., Kokkonen, P., Ledaine, I. and Niemkiewicz, E. 2006 Biovolumes and size-classes of phytoplankton in the Baltic Sea HELCOM Balt.Sea Environ. Proc. No. 106, 144pp <http://helcom.fi/Lists/Publications/BSEP106.pdf> (7 October 2015).

Utermöhl, H. 1958. Zur Vervollkommnung der quantitativen Phytoplankton-Methodik. *Mitt. int. Ver. theor. angew. Limnol.* 9: 1–38.



www.smhi.se



www.svenskalifewatch.se/en/



www.lifewatch.eu